Úlceras de Pé Diabético

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References: 1. INLOW's 60-second Diabetic Foot Screen. Screening tool. Canadian Association of Wound Care. www.cawc.net. 2011. 2. World Union of Wound Healing Societies (WUWHS). Florence Congress. Position Document. Local Management of Diabetic Foot Ulcers. Wounds International 2016. Available at: <http://www.woundsinternational.com> (Accessed 8 November 2018). 3. Frykberg, R.G., Banks, J. Challenges in the treatment of chronic wounds. Adv Wound Care (New Rochelle) 2015;4:560-82. 4. Ousey, K., Chadwick, P., Javien, A., et al. Identifying and treating foot ulcers in patients with diabetes: saving feet, legs and lives. J Wound Care 2018;27 (Suppl 5):S1-S52. 5. International Best Practice Guidelines. Wound Management in Diabetic Foot Ulcers. Wounds International 2013. Available at: <http://www.woundsinternational.com> (Accessed 8 November 2018). 6. White R. Evidence for atraumatic soft silicone wound dressing use. Wounds UK 2005;1(3):104-109. 7. Waring, M., Biefeldt, S., Matzold, K.P., Butcher, M. An evaluation of the skin stripping of wound dressing adhesives. J Wound Care 2011;20:412-22. 8. White, R. A multinational survey of the assessment of pain when removing dressings. Wounds UK 2008;4:14-22. 9. David, F., Wutze, J.-L., Bréton, N., et al. A randomized, controlled, non-inferiority trial comparing the performance of a soft silicone-coated wound contact layer (Mepitel One) with a lipidic oiled wound contact layer (UgoTul) in the treatment of acute wounds. Int Wound J 2017;doi:10.1111/iwj.12853. 10. Gotschall, C.S., Morrison, M.I., Eichelberger, M.R. Prospective, randomized study of the efficacy of Mepitel on children with partial-thickness scalds. J Burn Care Rehabil 1998;19:279-83. 11. Silverstein, P., Heimbach, D., Meites, H., et al. An open, parallel, randomized, comparative, multicenter study to evaluate the cost-effectiveness, performance, tolerance, and safety of a silver-containing soft silicone foam dressing (intervention) vs silver sulfadiazine cream. J Burn Care Res 2011;32:617-26. 12. Gee Kee, E.L., Kimble, R.M., Cuttle, L., Khan, A., Stockton, K.A. Randomized controlled trial of three burn dressings for partial thickness burns in children. Burns 2015;41:946-55. 13. Meaume, S., Van De Looverbosch, D., Heyman, H., Romanelli, M., Ciangherotti, A., Charpin, S. A study to compare a new self-adherent soft silicone dressing with a self-adherent polymer dressing in stage II pressure ulcers. Ostomy Wound Manage 2003; 49 (9): 44-51. 14. Wiberg, A.-B., Feili, F., Daun, E.-K. Preventing maceration with a soft silicone dressing: in-vitro evaluation. Poster presentation at the 3rd Congress of the World Union of Wound Healing Societies, Toronto, Canada, 2008. 15. Bredow, J., Hoffmann, K., Hellmich, M., Eysel, P., Zarghooni, K. Randomized clinical trial to evaluate performance of flexible self-adherent absorbent dressing coated with silicone layer after hip, knee or spinal surgery in comparison to standard wound dressing. Poster presentation at the 5th Congress of the World Union of Wound Healing Societies, Florence, Italy, 2016. 16. Patton, M.L., Mullins, R.F., Smith, D., Korentager, R. An open, prospective, randomized pilot investigation evaluating pain with the use of a soft silicone wound contact layer vs. traditional adhesive dressings. J Burn Care Res 2013;34:674-81. 17. Molnlycke Health Care data on file, report 20160105-002. 18. Chadwick, P., Taherinejad, F., Hamberg, K., Waring, M. Clinical and scientific data on a silver-containing soft-silicone foam dressing: an overview. J Wound Care 2009;18:483-91. 19. ProDerm study report 16.0180-23. Assessment of Wearing Properties of Wound Dressings on the Elbows. Data on file. 20. ProDerm study report 16.0180-23. Assessment of Wearing Properties of Wound Dressings on the Elbows. Data on file. 21. ALTEN Finite Element Modelling simulation. Laboratory report no. PD-530246. 22. Haycocks, S., Chadwick, P., Davies, P. Case series: Mepilex Border Comfort in the treatment of diabetic foot ulcers with exudate. Diabetic Foot Journal 2018;21:265-71. 23. External test lab report SMTL15/4863/2. 24. External lab report: NAMS A 09C 29253 01/09C 29253 02. 25. Kies C.L., Murrah, C.P., Smith, K., et al. Achieving and sustaining zero. Preventing surgical site infections after isolated coronary artery bypass with saphenous vein harvest through implementation of a staff-driven quality improvement process. Dimensions Crit Care Nurs 2015;34:265-72. 26. Sinet, S., Beele, H., Saine, L., Sny, E., Henricks, B. Open, non-comparative, multicentre post market clinical follow-up investigation to evaluate performance and safety on pressure ulcers when using a gelling fibre dressing as intended. Poster Presentation at European Pressure Ulcer Advisory Panel Conference, 2015, Ghent, Belgium. 27. Chadwick, P., McCauley, J. Open, non-comparative, multicentre post clinical study of the performance and safety of a gelling fibre wound dressing on diabetic foot ulcers. J Wound Care 2016;25:290-300. 28. Molnlycke Health Care. Laboratory report 20080616-007. Data on file. 29. Molnlycke Health Care. Laboratory report 20080616-006. Data on file. 30. Grocott Patricia Clinical investigation Mepilex Transfer, Clinical investigation of a silicone dressing in product development phase in the palliative management of patients with pressure sores and malignant wounds, study id MIN101 UK, London UK, 2000. 31. External lab report: NAMS A 11C_51788_01. 32. Meuleneire, F. Management of diabetic foot ulcers using dressings with SafetaC: a review of case studies. Wounds UK 2008;4:16-30. 33. Tickle, J., Fletcher, J. Mextra Superabsorbent made easy. Wounds UK 2012;8:1-4.

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© 2019 Molnlycke Health Care AB. Todos os direitos reservados. HQIM000894